

Love the everyday details in life

ReSound GN
.....



w.e.f 1st November 2023

ReSound OMNIA

Now Full Family of Organic Hearing




ReSound Premium RIE Hearing Aids


iSolate nanotech - All Weather Proof Hearing Aids
(Made for iPhone + Made for Android with ASHA + ReSound Assist Live)

Product Image	Model Name	No. of Channels	Receiver Options	Fitting Range	Warranty (Standard+ Promotional Yrs)	MRP Rs
---------------	------------	-----------------	------------------	---------------	--------------------------------------	--------


OMNIA MiniRIE Rechargeable Kit (M&RIE / Standard Receiver)

 <p>Li-ion</p>	OMNIA 960 MiniRIE (2 Hearing Aids + 1 Premium Charger)	17	MM/LP/MP/HP/UP	20-110	2+2	7,99,995
	OMNIA 960 MiniRIE (1 Hearing Aid + 1 Standard Charger)	17	MM/LP/MP/HP/UP	20-110	2+2	4,09,995
	OMNIA 760 MiniRIE (2 Hearing Aids + 1 Premium Charger)	14	MM/LP/MP/HP/UP	20-110	2+2	4,39,995
	OMNIA 760 MiniRIE (1 Hearing Aid + 1 Standard Charger)	14	MM/LP/MP/HP/UP	20-110	2+2	2,29,995
	OMNIA 560 MiniRIE (2 Hearing Aids + 1 Premium Charger)	12	MM/LP/MP/HP/UP	20-110	2+1	2,49,995
	OMNIA 560 MiniRIE (1 Hearing Aid + 1 Standard Charger)	12	MM/LP/MP/HP/UP	20-110	2+1	1,34,995
	OMNIA 460 MiniRIE (2 Hearing Aids + 1 Premium Charger)	12	MM/LP/MP/HP/UP	20-110	2	1,69,995
OMNIA 460 MiniRIE (1 Hearing Aid + 1 Standard Charger)	12	MM/LP/MP/HP/UP	20-110	2	94,995	

OMNIA Non-Rechargeable RIE (M&RIE / Standard Receiver)

 <p>Zinc Air 312</p>	OMNIA 961 RIE	17	MM/LP/MP/HP/UP	20-110	2+2	3,79,995
	OMNIA 761 RIE	14	MM/LP/MP/HP/UP	20-110	2+2	1,89,995
	OMNIA 561 RIE	12	MM/LP/MP/HP/UP	20-110	2+1	1,14,995
	OMNIA 461 RIE	12	MM/LP/MP/HP/UP	20-110	2	79,995

ONE 4 Rechargeable RIE Kit (Standard Receiver)

 <p>Li-ion</p>	ONE 461 RIE (2 Hearing Aids + 1 Premium charger)	12	LP/MP/HP/UP	20-110	2	1,19,995
	ONE 461 RIE (1 Hearing Aid + 1 Standard charger)	12	LP/MP/HP/UP	20-110	2	69,995

Warranty on chargers for rechargeable hearing aids is 1 year
Prices of Non - Rechargeable Hearing Aids are on per unit basis
:- 2.4 GHz direct wireless connectivity products Li-Ion stands for Lithium-Ion

Gain values as per IEC60118-0 Ear Simulator
Full on Gain at 50dB input: MM - 63, LP - 63, MP - 69, HP - 75, UP - 83

Improvement in speech understanding **150%**

Natural sound quality

90% of listeners prefer M&RIE


Don't worry about wind noise

15dB less wind noise

Easier phone calls

95% users reported effortless phone calls

Better ability to track moving sounds



ReSound Premium Custom Hearing Aids

iSolate nanotech - All Weather Proof Hearing Aids
(Made for iPhone + Made for Android with ASHA + ReSound Assist Live)

Product Pictures	Model Name	No. of Channels	Receiver Options	Fitting Range	Warranty (Standard+ Promotional Yrs)	MRP Rs
------------------	------------	-----------------	------------------	---------------	--------------------------------------	--------

OMNIA Rechargeable ITC Earbud Kit



OMNIA 9 ITC (2 Hearing Aids + 1 Custom Charger)		17	LP/MP/HP/UP	20-110	2+2	7,99,995
OMNIA 9 ITC (1 Hearing Aid + 1 Custom Charger)		17	LP/MP/HP/UP	20-110	2+2	4,09,995
OMNIA 7 ITC (2 Hearing Aids + 1 Custom Charger)		14	LP/MP/HP/UP	20-110	2+2	4,39,995
OMNIA 7 ITC (1 Hearing Aid + 1 Custom Charger)		14	LP/MP/HP/UP	20-110	2+2	2,29,995
OMNIA 5 ITC (2 Hearing Aids + 1 Custom Charger)		12	LP/MP/HP/UP	20-110	2+1	2,69,995
OMNIA 5 ITC (1 Hearing Aid + 1 Custom Charger)		12	LP/MP/HP/UP	20-110	2+1	1,44,995
OMNIA 4 ITC (2 Hearing Aids + 1 Custom Charger)		12	LP/MP/HP/UP	20-110	2	1,79,995
OMNIA 4 ITC (1 Hearing Aid + 1 Custom Charger)		12	LP/MP/HP/UP	20-110	2	99,995

ONE 4 Rechargeable ITC Earbud Kit

ONE 4 ITC (2 Hearing Aids + 1 Custom Charger)		12	LP/MP/HP/UP	20-110	2	1,41,995
ONE 4 ITC (1 Hearing Aid + 1 Custom Charger)		12	LP/MP/HP/UP	20-110	2	80,995

OMNIA Non Wireless CIC



OMNIA 9 CIC		17	LP/MP/HP	20-110	2+2	3,89,995
OMNIA 7 CIC		14	LP/MP/HP	20-110	2+2	2,09,995
OMNIA 5 CIC		12	LP/MP/HP	20-110	2+1	1,19,995
OMNIA 4 CIC		12	LP/MP/HP	20-110	2	74,995

ONE 4 Non Wireless CIC

ONE 4 CIC		12	LP/MP/HP	20-110	2	56,995
-----------	--	----	----------	--------	---	--------

OMNIA Non Wireless IIC (CIC with black face plate)



OMNIA 9 IIC		17	LP/MP	20-100	2+2	389,995
OMNIA 7 IIC		14	LP/MP	20-100	2+2	209,995
OMNIA 5 IIC		12	LP/MP	20-100	2+1	119,995
OMNIA 4 IIC		12	LP/MP	20-100	2	74,995

ONE Non Wireless IIC (CIC with black face plate)

ONE 4 IIC		12	LP/MP	20-100	2	56,995
-----------	--	----	-------	--------	---	--------

Warranty on chargers for rechargeable hearing aids is 1 year

Prices of Non - Rechargeable Hearing Aids are on per unit basis

:- 2.4 GHz direct wireless connectivity products Li-Ion stands for Lithium-Ion

Gain values as per IEC60118-0 Ear Simulator

ITC Rechargeable: Full on Gain at 50dB input: LP - 50, MP - 60, HP - 70, UP - 78

Gain values as per IEC60118-0 2cc Coupler


CIC Non-Rechargeable: Full on Gain at 50dB input: LP - 40, MP - 50, HP - 60

ReSound Premium Rechargeable BTE Hearing Aids

iSolate nanotech - All Weather Proof Hearing Aids
(Made for iPhone + Made for Android with ASHA + ReSound Assist Live)

Product Pictures	Model Name	No. of Channels	Receiver Options	Fitting Range	Warranty (Standard + Promotional Yrs)	MRP Rs
------------------	------------	-----------------	------------------	---------------	---------------------------------------	--------

OMNIA Rechargeable 88 BTE Kit




OMNIA 988 BTE (2 Hearing Aids + 1 Desktop Charger)	📶	17	HP BTE	30-110	2+2	7,99,995
OMNIA 988 BTE (1 Hearing Aid + 1 Desktop Charger)	📶	17	HP BTE	30-110	2+2	4,07,995
OMNIA 788 BTE (2 Hearing Aids + 1 Desktop Charger)	📶	14	HP BTE	30-110	2+2	4,39,995
OMNIA 788 BTE (1 Hearing Aid + 1 Desktop Charger)	📶	14	HP BTE	30-110	2+2	2,27,995
OMNIA 588 BTE (2 Hearing Aids + 1 Desktop Charger)	📶	12	HP BTE	30-110	2+1	2,49,995
OMNIA 588 BTE (1 Hearing Aid + 1 Desktop Charger)	📶	12	HP BTE	30-110	2+1	1,32,995
OMNIA 488 BTE (2 Hearing Aids + 1 Desktop Charger)	📶	12	HP BTE	30-110	2	1,41,995
OMNIA 488 BTE (1 Hearing Aid + 1 Desktop Charger)	📶	12	HP BTE	30-110	2	78,995

ONE 4 Rechargeable 88 BTE Kit

ONE 488 BTE (2 Hearing Aids + 1 Desktop Charger)	📶	12	HP BTE	30-110	2	1,13,995
ONE 488 BTE (1 Hearing Aid + 1 Desktop Charger)	📶	12	HP BTE	30-110	2	59,995

OMNIA Rechargeable 77 BTE Kit (Classic Hook/Thin Tube)



OMNIA 977 BTE (2 Hearing Aids + 1 Desktop charger)	📶	17	NP BTE	20-100	2+2	7,99,995
OMNIA 977 BTE (1 Hearing Aid + 1 Desktop charger)	📶	17	NP BTE	20-100	2+2	4,07,995
OMNIA 777 BTE (2 Hearing Aids + 1 Desktop charger)	📶	14	NP BTE	20-100	2+2	4,39,995
OMNIA 777 BTE (1 Hearing Aid + 1 Desktop charger)	📶	14	NP BTE	20-100	2+2	2,27,995
OMNIA 577 BTE (2 Hearing Aids + 1 Desktop charger)	📶	12	NP BTE	20-100	2+1	2,49,995
OMNIA 577 BTE (1 Hearing Aid + 1 Desktop charger)	📶	12	NP BTE	20-100	2+1	1,32,995
OMNIA 477 BTE (2 Hearing Aids + 1 Desktop charger)	📶	12	NP BTE	20-100	2	1,41,995
OMNIA 477 BTE (1 Hearing Aid + 1 Desktop charger)	📶	12	NP BTE	20-100	2	78,995

ONE 4 Rechargeable 77 BTE Kit (Classic Hook/Thin Tube)

ONE 477 BTE (2 Hearing Aids + 1 Desktop Charger)	📶	12	NP BTE	20-100	2	1,13,995
ONE 477 BTE (1 Hearing Aid + 1 Desktop Charger)	📶	12	NP BTE	20-100	2	59,995

Gain value as per IEC-60118-0-IEC-711-Ear Simulator

Full on Gain at 50dB input: NP BTE- 67, HP BTE- 80

Warranty on chargers for rechargeable hearing aids is 1 year

Prices of Non - Rechargeable Hearing Aids are on per unit basis

📶 :- 2.4 GHz direct wireless connectivity products

Li-Ion stands for Lithium-Ion



Know the Difference

	Speech Understanding	Noise Reduction	Automatic adjustments as per environment	Induction/Wireless charging	Apple & Android audio Streaming	Remote care/ Remote audiology care
OMNIA 9	★★★★★	★★★★★	★★★★★			
OMNIA 7	★★★★	★★★★	★★★★			
OMNIA 5	★★★	★★★	★★★			
OMNIA 4	★★	★★	★★			
ONE 4	★	★	★			

Speech Understanding

Feature	Benefit	OMNIA 9	OMNIA 7	OMNIA 5	OMNIA 4	ONE 4
WARP compression channels	Replicates the way the natural ear process sound "Near-zero distortion & remarkable purity of sound"	17	14	12	12	12
All-Around 360	Hearing clarity from all directions, provides brain with natural acoustic information "Industry best Environment mapping"	✓	--	--	--	--
Front Focus	Hear best in noise "Upto 150% better speech understanding in noise"	✓	--	--	--	--
All Access Directionality	Automatically tune into what you like to hear with "Artificial Intelligence"	✓	✓	--	--	--
Ultra Focus	Easy focus on speech coming from front (almost 4.3 dB more focus) "Adaptive Beamforming"	✓	✓	--	--	--
Binaural Directionality III	Hear both environment sounds & speech, without compromising each other "Automatic Balanced Hearing"	✓	✓	✓	--	--
Spatial Sense	Locate the sounds precisely front vs back & Up vs Down with "Inter ear Compression" & "Pinna Restoration"	✓	✓	✓	--	--
Binaural Directionality	Reduced hearing effort with "Industry's fastest Device to Device communication"	✓	✓	✓	✓	--
M&RIE	Uses unique ear shape to deliver natural sound, easily locate sounds & wind noise reduction. "Industry's first & Only 3 Mic technology"	✓	✓	✓	✓	--
Natural Own Voice Perception	Own Voice is processed differently from others voice, reducing over amplification of own voice "Band Split processor"	✓	✓	✓	✓	--
Synchronised Soft Switching	Automatic focus on speech in noisy situations "Synchronous auto-adaptation" of hearing aids in both ears	✓	✓	✓	✓	✓
Sound Shaper	Hear inaudible sounds back again even if you have high hearing loss/ dead regions "Linear Frequency Compression"	✓	✓	✓	✓	✓

Noise Reduction

Noise Tracker II	Reduce unwanted sounds without degrading speech "Spectral subtraction Algorithms"	10dB Noise Reduction	8dB Noise Reduction	6dB Noise Reduction	3dB Noise Reduction	3dB Noise Reduction
Wind Guard	Reduces disturbing wind noise sounds without reducing speech "Twin-mic correlation coefficient Algorithms"	Strong	Moderate	Mild	Mild	Mild
Personalized Noise Cancellation	Reduce noise of only specific sounds of an environment that you don't like "Per environment Noise reduction"	✓	--	--	--	--
Environmental Optimiser	Confidently change environments, your hearing aids will adapt automatically "Industry's best 98% accurate environmental classifier"	✓	✓	✓	--	--
DFS Ultra	Rich sound without annoying whistling sounds "Feedback path modeling" blocks whistling sounds	Strong	Strong	Strong	Strong	Strong

Comfort & Convenience


Synchronized Acceptance Manager	Intelligent technology to get you accustomed to new hearing aids "Data-Logging technology" ensures, your hearing aids track your progress	✓	✓	✓	✓	✓
Tinnitus Sound Generator	Hearing & Ringing sound in the ear have common origin & can be treated with TSG built using "Sound & Cognitive behaviour therapy" principles	✓	✓	✓	✓	✓
Ear to Ear Communication	Control both of your hearing aids with a single touch "Industry unique 2.4Ghz Ear-to-Ear technology"	✓	✓	✓	✓	✓

Range of World's Best Super Power Hearing Aid's just got Bigger

ReSound ENZO Q (World's Best Super Power Hearing Aid)


iSolate nanotech - All Weather Proof Hearing Aids with 2 Years Warranty

(Made for iPhone + Made for Android with ASHA + ReSound Assist Live)


Product Pictures	Model Name	No. of Channels	Full on Gain 50dB	Fitting Range	Warranty	MRP Rs
	ENZO Q 998 SP BTE	17	86	60-120	2	3,49,995
	ENZO Q 798 SP BTE	14	86	60-120	2	1,49,995
	ENZO Q 598 SP BTE	12	86	60-120	2	69,995

Buy 2 ReSound ENZO Quattro Hearing Aids of MRP 1,49,995/- each or higher & get Any 1 Wireless Accessory worth MRP 17,700/- FREE


ReSound OMNIA (World's Most Powerful Rechargeable PBTE)

Product Pictures	Model Name	No. of Channels	Full on Gain 50dB	Fitting Range	Warranty (Standard+ Promotional Yrs)	MRP Rs
OMNIA Rechargeable 88 BTE Kit						
	OMNIA 988 BTE (2 Hearing Aids + 1 Desktop Charger)	17	80	30-110	2+2	7,99,995
	OMNIA 988 BTE (1 Hearing Aid + 1 Desktop Charger)	17	80	30-110	2+2	4,07,995
	OMNIA 788 BTE (2 Hearing Aids + 1 Desktop Charger)	14	80	30-110	2+2	4,39,995
	OMNIA 788 BTE (1 Hearing Aid + 1 Desktop Charger)	14	80	30-110	2+2	2,27,995
	OMNIA 588 BTE (2 Hearing Aids + 1 Desktop Charger)	12	80	30-110	2+1	2,49,995
	OMNIA 588 BTE (1 Hearing Aid + 1 Desktop Charger)	12	80	30-110	2+1	1,32,995
	OMNIA 488 BTE (2 Hearing Aids + 1 Desktop Charger)	12	80	30-110	2	1,41,995
	OMNIA 488 BTE (1 Hearing Aid + 1 Desktop Charger)	12	80	30-110	2	78,995
ONE 4 Rechargeable 88 BTE Kit						
	ONE 488 BTE (2 Hearing Aids + 1 Desktop Charger)	12	80	30-110	2	1,13,995
	ONE 488 BTE (1 Hearing Aid + 1 Desktop Charger)	12	80	30-110	2	59,995

ReSound Key (World's Only 2.4GHz Wireless SP BTE)

Product Pictures	Model Name	No. of Channels	Full on Gain 50dB	Fitting Range	Warranty	MRP Rs
	Key 498 SP BTE	12	86	60-120	2	48,995
	Key 398 SP BTE	8	86	60-120	2	37,995
	Key 298 SP BTE	6	86	60-120	2	27,995

ReSound Magna

Product Pictures	Model Name	No. of Channels	Full on Gain 50dB	Fitting Range	Warranty	MRP Rs
	Magna 490 SP BTE	9	86	60-120	2	22,495

Warranty on chargers for rechargeable hearing aids is 1 year
Prices of Non - Rechargeable Hearing Aids are on per unit basis

:- 2.4 GHz direct wireless connectivity products Li-Ion stands for Lithium-Ion Gain value as per IEC-60118-O-IEC-711- Ear Simulator



ReSound OMNIA HP BTE CHARGER

Know the Difference

	Speech Understanding	Noise Reduction	Automatic adjustments as per environment	Induction/Wireless charging	Apple & Android audio Streaming	Remote care/ Remote audiology care
ENZO 9	★★★★★	★★★★★	★★★★★			
ENZO 7	★★★★	★★★★	★★★★	★★★★★	★★★★★	★★★★★
ENZO 5	★★★	★★★	★★★			

Speech Understanding

Feature	Benefit	ENZO 9	ENZO 7	ENZO 5
WARP compression channels	Replicates the way the natural ear process sound "Near-zero distortion & remarkable purity of sound"	17	14	12
Binaural Directionality III	Hear both environment sounds & speech, without compromising each other "Automatic Balanced Hearing"	✓	--	--
Spatial Sense	Locate the sounds precisely front vs back & Up vs Down with "Inter ear Compression" & "Pinna Restoration"	✓	--	--
Autoscope Adaptive Directionality	Advanced form of Directionality automatically helps focus on speech in noisy situation "Focus on a person or group of people as needed"	✓	--	--
Binaural Directionality	Reduced hearing effort with "Industry fastest Device to Device communication"	✓	✓	--
Synchronised Soft Switching	Automatic focus on speech in noisy situations "Synchronus auto-adaptation" of hearing aids in both ears	✓	✓	✓
Sound Shaper	Hear inaudible sounds back again even if you have high hearing loss/ dead regions "Linear Frequency Compression"	✓	✓	✓

Noise Reduction

Noise Tracker II	Reduce unwanted sounds with out degrading speech "Spectral subtraction Algorithms"	10dB Noise Reduction	8dB Noise Reduction	6dB Noise Reduction
Wind Guard	Reduces disturbing wind noise sounds without reducing speech "Twin-mic correlation coefficient Algorithms"	Strong	Moderate	Mild
Personalized Noise Cancellation	Reduce noise of only specific sounds of an environment that you don't like "Per environment Noise reduction"	✓	--	--
Environmental Optimiser	Confidently change environments, your hearing aids will adapt automatically "Industry best 98% accurate environmental classifier"	✓	✓	--
DFS Ultra	Rich sound without annoying whistling sounds "Feedback path modeling" blocks whistling sounds	Strong	Strong	Strong

Comfort & Convenience

Synchronized Acceptance Manager	Intelligent technology to get you accustomed to new hearing aids "Data-Logging technology" ensures, your hearing aids track your progress	✓	✓	
Tinnitus Sound Generator	Hearing & Ringing sound in the ear have common origin & can be treated with TSG built using "Sound & Cognitive behaviour therapy" principles	✓	✓	✓
Ear to Ear Communication	Control both of your hearing aids with a single touch "Industry unique 2.4Ghz Ear-to-Ear technology"	✓	✓	✓

ReSound Key

iSolate nanotech - All Weather Proof Hearing Aids with 2 Years Warranty

(Made for iPhone + Made for Android with ASHA + ReSound Assist Live)

Product Pictures	Model Name	No. of Channels	Receiver Options	Fitting Range	MRP Rs
------------------	------------	-----------------	------------------	---------------	--------



ReSound Key Rechargeable RIE

(Made for iPhone + Made for Android with ASHA + ReSound Assist Live)

Key 461 RIE Kit (2 Hearing Aids + 1 Desktop Charger)	📶	12	LP/MP/HP/UP	20-110	99,995 NEW
Key 461 RIE	📶	12	LP/MP/HP/UP	20-110	49,995 NEW
Key 361 RIE Kit (2 Hearing Aid + 1 Desktop Charger)	📶	8	LP/MP/HP/UP	20-110	83,995 NEW
Key 361 RIE	📶	8	LP/MP/HP/UP	20-110	41,995 NEW



ReSound Key Non-Rechargeable RIE

Made for iPhone + Made for Android with ASHA + ReSound Assist Live (MFI & ASHA streaming only in Key 3 & 4)

Key 361 RIE	📶	8	LP/MP/HP/UP	20-110	36,995
Key 362 RIE	📶	8	LP/MP/HP/UP	20-110	36,995
Key 261 RIE	📶	6	LP/MP/HP/UP	20-110	27,495
Key 262 RIE	📶	6	LP/MP/HP/UP	20-110	27,495



ReSound Key Super Power BTE

Made for iPhone + Made for Android with ASHA + ReSound Assist Live (MFI & ASHA streaming only in Key 3 & 4)

Key 498 SP BTE	📶	12	SP BTE	60-120	48,995
Key 398 SP BTE	📶	8	SP BTE	60-120	37,995
Key 298 SP BTE	📶	6	SP BTE	60-120	27,995



ReSound Key Normal & High Power BTE

Made for iPhone + Made for Android with ASHA + ReSound Assist Live (MFI & ASHA streaming only in Key 3 & 4)


Key 488 HP BTE	📶	12	HP BTE	30-110	48,995
Key 477 NP BTE	📶	12	NP BTE	20-100	47,995
Key 388HP BTE	📶	8	HP BTE	30-110	37,995
Key 377 NP BTE	📶	8	NP BTE	20-100	36,995
Key 288 HP BTE	📶	6	HP BTE	30-110	28,495
Key 277 NP BTE	📶	6	NP BTE	20-100	27,995
Key 188 HP BTE	📶	4	HP BTE	30-110	20,995
Key 177 NP BTE	📶	4	NP BTE	20-100	20,495



ReSound Key Wireless ITC

Made for iPhone + Made for Android with ASHA + ReSound Assist Live (MFI & ASHA streaming only in Key 3 & 4)

Key 4 ITC	📶	12	MP/HP/UP	20-110	49,995
Key 3 ITC	📶	8	MP/HP/UP	20-110	38,995
Key 2 ITC	📶	6	MP/HP/UP	20-110	28,995



ReSound Key Wireless CIC

Made for iPhone + Made for Android with ASHA + ReSound Assist Live

Key 4 CIC	📶	12	LP/MP/HP	20-105	40,995
Key 3 CIC	📶	8	LP/MP/HP	20-105	29,995

Gain values as per IEC60118-0 Ear Simulator

RIE - Full on Gain at 50dB input: LP - 62, MP - 67, HP - 74, UP - 82

ITC/CIC - Full on Gain at 50dB input: LP - 49, MP - 59, HP - 69, UP - 76

BTE - Full on Gain at 50dB input: NP BTE - 66, HP BTE - 80

Warranty on chargers for rechargeable hearing aids is 1 year

Prices of Non - Rechargeable Hearing Aids are on per unit basis

📶 :- 2.4 GHz direct wireless connectivity products Li-Ion stands for Lithium-Ion



Desktop Charger

Know the Difference

	Speech Understanding	Noise Reduction	Automatic adjustments as per environment	Induction/Wireless charging	Apple & Android audio Streaming	Remote care/ Remote audiology care
KEY 4	★★★★★	★★★★★	★★★★★	★★★★★	★★★★★	★★★★★
KEY 3	★★★★★	★★★★★	★★★★★	★★★★★	★★★★★	
KEY 2	★★★	★★★	--	--	--	
KEY 1	★★	★★	--	--	--	

Speech Understanding



Feature	Benefit	KEY 4	KEY 3	KEY 2	KEY 1
WARP compression channels	Replicates the way the natural ear process sound "Near-zero distortion & remarkable purity of sound"	12	8	6	4
Natural Directionality II	User brain gets binaural image for better representation of the speech.doesn't require user to change programs. "Increases ease of listening"	✓	--	--	--
Synchronised Soft Switching	Automatic focus on speech in noisy situations "Synchron auto-adpation" of hearing aids in both ears	✓	--	--	--
Soft Switching	Smoothly switches the hearing aid between omni and directional mode depending on the sound environment. "Improves speech intelligibility in noisy situations"	✓	✓	--	--
Multiscope Directionality	Advanced form of Directionality that allows user to decide on how many people he/she wants to focus in noisy situation "Focus on a person or group of people as needed"	✓	✓	--	--
Adaptive Directionality	Dynamically changes the directional pattern to reduce the loudest sounds from the rear to. "Better concentrate on signal from the front"	✓	✓	✓	✓
Environmental Classifier	Analyzes the incoming sounds and classifies them into listening environments sets foundation for a "comfortable hearing experience"	✓	✓	✓	✓
Impulse Noise Reduction	Improves listening comfort and sound quality for impulse sounds like clanking silverware or jangling key "Loud sounds wont bother you"	✓	--	--	--
Noise Tracker II	Improves listening comfort "without compromising speech understanding"	✓	✓	✓	✓
Wind Guard	Accurately detects wind noise & adaptively reduces it with out reducing volume "Comfort in windy situaions"	✓	✓	✓	--
DFS Ultra II	Reduces Whistling sound,preserving amplification& sound quality "No whistling sound or sudden reduction of volume"	✓	✓	✓	✓
Expansion	Reduces amplification for very soft sounds which are not of interest to the wearer, such as refrigerator, electronic equipment humming etc. "Unwanted soft noise reduction"	✓	✓	✓	✓
Synchronized Acceptance Manager	Intelligent technology to get you accoutmed to new hearing aids. "Data-Logging technology" ensures, your hearing aids track your progress	✓	✓	--	--
Direct audio streaming (Made for i-Phone & Android ASHA)	High definition audio streaming/routing from your phones to your hearing aids. "Your hearing aids can work as bluetooth headphones"	✓	✓	--	--
ReSound Assist Live	Tele-audiology – Provide's individualised hearing care, in circumstances when patient cant reach clinic for adjustments	✓	✓	✓	✓
ReSound Smart 3D App	"Personalise and control your hearing aids" with just a swipe on your phone	✓	✓	✓	✓

Noise Reduction


Comfort & Convenience

ReSound LiNX Quattro Custom Hearing Aids



iSolate nanotech - All Weather Proof Hearing Aids with 2 Years Warranty
Made for i-Phone, Made for Andriod with ASHA, Resound Live Assist

Product Pictures	Model Name	No. of Channels	Receiver Options	Fitting Range	MRP Rs
	LiNX Quattro Wireless ITC				
	LiNX Q 9 ITC	17	LP/MP/HP/UP	20-110	3,24,995
	LiNX Q 7 ITC	14	LP/MP/HP/UP	20-110	1,69,995
	LiNX Q 5 ITC	12	LP/MP/HP/UP	20-110	1,09,995
	LiNX Quattro Wireless CIC				
	LiNX Q 9 CIC	17	LP/MP/HP/UP	20-110	3,24,995
	LiNX Q 7 CIC	14	LP/MP/HP/UP	20-110	1,69,995
	LiNX Q 5 CIC	12	LP/MP/HP/UP	20-110	1,09,995

ReSound LiNX 3D Custom Hearing Aids

Product Pictures	Model Name	No. of Channels	Receiver Options	Fitting Range	MRP Rs
	LiNX 3D Non Wireless CIC with UP receiver				
	LiNX 3D 9 CIC	17	UP	30-110	389,995
	LiNX 3D 7 CIC	14	UP	30-110	209,995
	LiNX 3D 5 CIC	12	UP	30-110	89,995

ReSound Enya Custom Hearing Aids

Product Pictures	Model Name	No. of Channels	Receiver Options	Fitting Range	MRP Rs
	Enya Non Wireless CIC				
	Enya 4 CIC	10	LP/MP/HP/UP	20-110	40,995
	Enya 3 CIC	8	LP/MP/HP/UP	20-110	33,300
	Enya 2 CIC	6	LP/MP/HP/UP	20-110	25,600
	Enya IIC (CIC with black face plate)				
	Enya 4 IIC	10	LP/MP	20-100	40,995
	Enya 3 IIC	8	LP/MP	20-100	33,300
	Enya 2 IIC	6	LP/MP	20-100	25,600




Gain values as per IEC60118-0 Ear Simulator

ITC/CIC/IIC - Full on Gain at 50dB input: LP - 49, MP - 59, HP - 69, UP - 79

📶 :- 2.4 GHz direct wireless connectivity products

ReSound Vea

iSolate nanotech - All Weather Proof Hearing Aids with 2 Years Warranty

Product Pictures	Model Name	No. of Channels	Receiver Options	Fitting Range	MRP Rs
VEA Non Wireless In The Canal Hearing Aids (ITC)					
 <p>Zinc Air 312</p>	Vea 3 ITC	7	S/P	15-90	22,600
	Vea 2 ITC	6	S/P	15-90	20,500
	Vea 1 ITC	4	S/P	15-90	16,495
VEA Non Wireless Completely In the Canal Hearing Aids (CIC)					
 <p>Zinc Air 10 A</p>	Vea 3 CIC	7	S/P	10-90	22,600
	Vea 2 CIC	6	S/P	10-90	20,500
	Vea 1 CIC	4	S/P	10-90	16,495
VEA Normal BTE					
 <p>Zinc Air 13</p>	Vea 170 VI BTE <small>Thin tube/Closed</small>	4	NP	10-110	14,495

Gain values as per IEC60118-0 Ear Simulator
 ITC/CIC - Full on Gain at 50dB input: S - 52, P - 60
 BTE - Full on Gain at 50dB input: NP - 69

ReSound Vea Sounds Good



Accessories MRP Inclusive of GST

Sr. No.	Chargers for Rechargeable Hearing Aids (with 1 year warranty)	MRP Per Unit (Rs)
1	RIE Premium Charger (ReSound OMNIA, ONE, Key)	15,495
2	RIE Standard Charger (ReSound OMNIA, ONE)	10,495
3	Custom Charger (ReSound OMNIA ITC, ONE ITC)	10,495
4	Desktop Charger (ReSound OMNIA BTE, ONE BTE, Key RIE)	8,495

Sr. No.	Wireless Accessories (with 1 year warranty)	MRP Per Unit (Rs)
5	Phone Clip+	17,700
6	TV Streamer 2	17,700
7	Remote control 2	17,700
8	Multi Mic	17,700

Sr. No.	Accessories		MRP Per Unit (Rs)
9	M&RIE Receivers	Available in 0, 1, 2, 3 & 4 Sizes	9,900
10	RIE Receivers (excluding M&RIE receivers)	Available in 0, 1, 2, 3 & 4 Sizes	5,900
11	M&RIE Microphone Filters	Pack includes 8 filter pieces	700
12	Rechargeable BTE Microphone Filter	Pack includes 8 filter pieces	700
13	Micromould For RIE Products		1,400
14	Hard BTE Earmould		1,500
15	Soft BTE Earmould		1,400
16	HF3 Filter	Strip include 15 filter pieces	1,050
17	Cerustop Waxguard Box	Strip include 8 filter pieces	700
18	Reshelling Charges	Change in canal size, damaged etc.	4,000
19	ReSound zinc air batteries (Pack of 6)	Battery size : 675,13,312,10A	330

Sr. No.	Extended Warranty		MRP Per Unit (Rs)
19	1 Year Extended Warranty	Applicable on MRP > Rs.15,000/-	10% of MRP + 18% GST
20	2 Year Extended Warranty	Applicable on MRP > Rs.15,000/-	15% of MRP + 18% GST



Multi Mic



TV Streamer 2



Remote Control 2



Phone Clip +

Charging options for ReSound hearing aids

ReSound offers rechargeable solutions with inductive charging and long battery life to keep your clients on the go.



Premium charger

Standard charger

Desktop charger


Custom charger

ReSound OMNIA™ miniRIE		06	06	-	-
ReSound OMNIA™ RIE ReSound ONE™ RIE		01	01	01	-
ReSound OMNIA™ BTE ReSound ONE™ BTE		-	-	70	-
ReSound OMNIA™ PBTE ReSound ONE™ PBTE		-	-	80	-
ReSound OMNIA™ ITE/ITC ReSound ONE™ ITE/ITC		-	-	-	C
ReSound LiNX Quattro™ RIE ReSound Key™ RIE		03	-	03	-

How to match ReSound hearing aids with the correct charger


Each charger type is offered with different inserts – known as the charging form factor (CFF) – which relate to a specific ReSound hearing aid collection or style.

The CFF is indicated with a number and is shown inside the charger and on the charger packaging.



Easier onboarding and better wearing experience with “Check my fit”

70% of users wore their hearing aids better with ‘Check My Fit’ improving sound experience by 11dB





Style meets Hearing

Rechargeable
Earbud

ReSound promises shell replacement
free for lifetime*

*Terms and conditions apply






Sounds natural

Hear in noise and all around. Now you can hear all the details you couldn't before, especially in noisy situations.



Feels natural

Forget you're wearing hearing aids. Our best wearing comfort in a discreet, comfortable design.



Connects you naturally

Stream audio, calls and more. Stay connected with greater streaming straight to your ears without interruptions.

ReSound OMNIA



Introducing world's smallest rechargeable RIE (MiniRIE)

GN Hearing India Private Limited
NMS Titanium, 4th Floor, Office No.401- 408, Plot No. 74, Sector 15, CBD Belapur, Navi Mumbai - 400614
Customer Helpline No. : +91 22 2752 4700 | www.resound.com/en-IN | www.pro.resound.com/en-in

www.facebook.com/resoundhearingindia | www.linkedin.com/company/resound

© 2021 GN Hearing Care Corporation. All rights reserved. ReSound is a trademark of GN Hearing Care Corporation. Apple, the Apple logo, iPhone, iPad, iPod touch and Apple Watch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries. Android, Google Pixel, Google Play and the Google Play logo are trademarks of Google Inc. Dolby and the double-D symbol are registered trademarks of Dolby Laboratories.

*Only when 2 hearing aids are worn

ReSound GN